

DUKE EXTRA 2.0 ACTIVE MATTRESS



Product Code: HAMDUKE002 (1200mm wide)

Description:

- Built with the same technologies as the Duke replacement mattress system. This mattress has been designed to fit the wider beds in our product range such as the Community and the Ultimate 1200 wide Profiling Beds.
- 29 cells (25 alternating cells with 4 static cells for the pillow section)
- With a maximum patient weight of 50 Stone, this mattress is suitable for larger patients

Mattress Dimensions: 2000 x 1200 x 200mm

Weight: 20kg

Max. Patient Weight: Harvest Pump 3: 50 Stone (317 kg)

Pressure Relief Risk: Very High

Pump Technical Specification: Please see separate pump specification sheet

Cleaning Guidance:

Cover should be wiped down with hot water and mild detergent. If heavily soiled, the cover may be machine washed up to 85°C and tumble dried. For disinfection use sodium hypochlorite solution max 1000ppm. Use no bleach, biological or phenolic-based cleaners. Dry completely before storing.

Warranty

In terms of Harvest Healthcare Guarantee and Warranties the Harvest Duke Extra 2.0 mattress is guaranteed for three years from date of purchase.

Exclusions:

1. Removable covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded
3. Components of pumps which are recommended for annual service, are guaranteed for one year.

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact.
(Mattress risk categories are to be used as a guide only)

Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 1: Feb.2020